ENGLISH

Intended use	For in vitro diagnostic use.
	Monoclonal Mouse Anti-Cytomegalovirus, Clones CCH2 + DDG9, is intended for use in immunohistochemistry. The antibody labels human cytomegalovirus (HCMV)-infected cells and is a useful aid in the identification of HCMV infection in human tissues (1, 2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Human cytomegalovirus is a ubiquitous herpes virus. The primary infection usually occurs in childhood where it most often is asymptomatic. After initial exposure, the virus establishes a life-long infection. In neonates, and in immunocompromised individuals, such as AIDS patients and transplant recipients, HCMV can cause severe life-threatening infections, and it is the most common cause of viral birth defects (5, 6). The virus retains the capability of reactivation which may occur following blood transfusions, pregnancy, solid-organ or bone marrow transplantation, immunosuppressive therapy, or other viral infections (3, 5). HCMV infects leukocytes, endothelial cells, connective tissue cells and epithelial cells, and is transmitted through most bodily fluids (5). Members of the herpes virus family have been implicated in the etiology of several human cancers. In the case of HCMV it has been associated with cervical carcinoma, adenocarcinomas of the prostate and colon, and Kaposi's sarcoma. HCMV infections have been shown to modulate the expression of various proteins involved in cell cycle regulation and apoptosis (5).
Reagent provided	A mixture of two monoclonal mouse antibodies provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . Clones: CCH2 and DDG9 (3). Isotypes: IgG1, kappa (CCH2), and IgG2a, kappa (DDG9). Mouse IgG concentration: see label on vial.
Immunogen	Crude cell lysates from cytomegalovirus strain Ad169-infected human embryonic fibroblasts (3).
Specificity	Antibody CCH2 reacts with an early nuclear protein identical with the non-structural DNA-binding protein p52 (7). More than 200 human CMV isolates have been successfully typed by this antibody which, therefore, seems to recognize an epitope which is highly conserved among HCMV strains (3, 4). Antibody DDG9 reacts with an immediate early nuclear protein of about 76 kDa (3). For both antibodies, the reactivity persists also at later stages during the HCMV infection where the localisation is less distinctly nuclear and appears to be in the cytoplasm. Laser confocal microscopy, however, shows that the reaction is limited to the nuclear membrane (3). The antibodies do not cross-react with adenovirus, herpes simplex virus, and varicella zoster virus (3).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. This product contains material of animal origin, therefore the product should be handled as potentially infectious. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue or cell pellet sections fixed in formalin. Pre-treatment of tissues with trypsin (1) or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval optimal results are obtained with Dako Target Retrieval Solution, High pH, Code S3308. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and pre-treatment of tissues with proteinase K. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. Frozen sections and cell preparations: The antibody can be used for labelling infected cell preparations fixed for 10 minutes in cold methanol (3, 4).
Staining procedure	Dilution: Monoclonal Mouse Anti-Cytomegalovirus, Code M0854, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of HCMV-infected human dermal fibroblasts and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, High pH, Code S3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is a mixture of (2 parts Dako Mouse IgG1, Code X0931, 1 part Dako Mouse IgG2a, Code X0943, and 1 part diluent) diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. Positive and negative controls should be run simultaneously with patient specimen. Visualization: Dako LSAB™+/HRP kit, Code K0679, and Dako EnVision™+/HRP kits, Codes K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, Code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.
Performance characteristics	Cells labelled by the antibody display a nuclear staining pattern early in the HCMV infection. At a later stage, an apparent cytoplasmic staining pattern might be observed (3).

Normal tissues: Screening of normal human tissues from donors without HCMV infection with the CCH2 antibody showed no reaction with heart, lung, liver, kidney, spleen, lymph node, tonsil, bone marrow, skeletal muscle, brain, thyroid gland, parathyroid gland, pancreas and adrenal gland (1).

Abnormal tissues: In HCMV-infected tissues, the CCH2 antibody labelled cells in lung, liver, kidney, lymph node, brain, thyroid and parathyroid glands (1).

The CCH2 and DDG9 antibodies together labelled HCMV-infected tissue of oesophagus, stomach, small intestine, colon, gall bladder, lung, adrenal gland, ovary and neural tissue (2).

The CCH2 antibody demonstrated a sensitivity for detection of HCMV infection of 85% and a specificity of 96% when applied on cultured cells 24 hours after inoculation with urine samples from 160 kidney transplanted patients. The sensitivity was comparable to that obtained by virus isolation (4). The sensitivity for detection of HCMV early antigen in cell culture after short time incubation increases when Antibody CCH2 and Antibody DDG9 are used in combination (3).

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1. Niedobitek G, Finn T, Herbst H, Gerdes J, Grillner L, Landqvist M, et al. Detection of cytomegalovirus by in situ hybridisation and immunohistochemistry using new monoclonal antibody CCH2: a comparison of methods. J Clin Pathol 1988;41:1005-9.
2. Saiz E, Lubin J, Robinson MJ. The modified Steiner stain: a new use for an old stain? Staining cytomegalovirus-infected cells in gastrointestinal biopsies. Histochem J 1998;30:549-52.
3. Wirgart BZ. Development of rapid methods for early diagnosis of cytomegalovirus infections [dissertation]. Stockholm:Karolinska Hospital; 1991.
4. Wirgart BZ, Landqvist M, Hökeberg I, Eriksson B-M, Olding-Stenkist E, Grillner L. Early detection of cytomegalovirus in cell culture by a new monoclonal antibody, CCH2. J Virol Methods 1990;27:211-20.
5. Doniger J, Muralidhar S, Rosenthal LJ. Human cytomegalovirus and human herpesvirus 6 genes that transform and transactivate. Clin Microbiol Rev 1999;12:367-82.
6. Jarvis MA, Nelson JA. Human cytomegalovirus persistence and latency in endothelial cells and macrophages. Curr Opin Microbiol 2002;5:403-7.
7. Plachter B, Nordin M, Wirgart BZ, Mach M, Stein H, Grillner L, et al. The DNA-binding protein P52 of human cytomegalovirus reacts with monoclonal antibody CCH2 and associates with the nuclear membrane at late times after infection. Virus Res 1992;24:265-76.

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